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Function and indication

All IST® appliances from Prof. Hinz have **two functions**:

1. To fix the lower jaw during sleep in a protruded position and
2. at the same time, reposition the tongue, tongue base and adjacent soft tissue holding all in forward position.

This helps to:

- achieve an enlargement of the rear pharyngeal space
- prevent the rear airway from being occluded or obstructed by a collapsing tongue.

The IST® appliance is designed to:

- reduce or avoid primary – harmless – snoring, and
- prevent temporary respiratory arrest in cases of a mild obstructive sleep apnoea syndrome.

Advantages of the IST® appliance:

- thin thermoformed splints do not restrict tongue movement
- unlimited mouth opening movements
- undisturbed mouth breathing
- lateral movements are possible
- good retention without wire elements
- individual adjustment of the lower jaw forward displacement by the dentist

1. Occlusion of rear airway by the collapsing tongue.



2. Mechanical forward displacement of the lower jaw by a IST® appliance including the tongue.

Preparations

1. In order to obtain the master models, these should be duplicated or the existing one impression should be cast 2x.
2. The working models shall be articulated with the aid of the IST® bite fork in the therapeutic forward position. With the help of the parallelometer, the prosthetic equator of each tooth is determined and marked.
3. Undertrained sites such as pontics, switching gaps or interdental spaces are to be blocked out with wax.
4. Before all thermoforming operations, when embedding the models in the granules, care should be taken that a distance of about 3 to 4 mm from the sulcus to the granulate is maintained.
5. The insulating film is deep drawn according to the data. In advance of the subsequent deep drawing, it is advisable to cut the insulating film interdentally by means of a scalpel and trim the film around the model.
6. A 2 mm thick film is now pulled over the prepared model e.g. either Duran or LAMlone.

The IST[®] classic appliance

- 1 Deep draw splints for the upper and lower jaw
- 2 Telescopes
- 3 Incorporated fixing elements
- 4 Protrusion screw for infinitely variable adjustment
- 5 Lateral bite planes
- 6 Screw for fixing of telescopes

Fitting of anchorage parts


1. Determination of upper/lower anchorage points
2. Fixing of anchorage parts with impression silicone
3. Fitting of anchorage parts with acrylic


Applying bite plane with acrylic in the upper and lower jaw

Under no circumstances lift the foils from the model!

- Finishing acrylic on model

Checking the occlusion in articulator with approx. 2 mm distance between the acrylic bite planes and more than 2 mm in the molar area.

Pressure forming of LAMItwo film

(Only with SCHEU appliances: Biostar®)

- Embedding model in granulate
- The model should be embedded as deep as possible (see fig.)
- Clamp the LAMItwo film with the soft side up
- Conditioning of the film with the enclosed bonding agent cloth. Subsequently, the film is deep-drawn

Finishing and final mounting*

- Lifting off splints from models and finishing
- Fitting of metal parts


* The necessary working tools for finishing are offered in a set (Art. no. 98 700).

The completed IST®classic appliance

Clean the IST®classic appliance on a daily basis using a toothbrush.



The IST[®] classic^{NEW} appliance

- 1 Deep drawing splints
- 2 IST[®]classic^{NEW} guiding bridge
- 3 Threaded rod (infinitely variable adjustment)
- 4 Lateral bite plane in the upper and lower jaw

ATTENTION: New fitting position!

Adjustment of the anchorage elements

1. After adjusting the protrusion guiding rods with acrylate, the anchorage elements are attached to the upper jaw splint and polymerised in the pressure pot.
2. The anchorage elements are adjusted and fixed onto the LAMlone film in the molar area of the lower jaw.
3. After adjusting the protrusion guiding rods, the anchorage elements are fixed with acrylic on the lower jaw splint and polymerised in the pressure pot.
4. Side view of the adjusted anchorage elements in the articulator.
5. Adjustment of the protrusion guiding rods takes place in the articulator.
6. **To maximise protrusion, the thread rod has to be shortened.**

Lamination technique with LAMItwo

1

1. The LAMItwo film is put in the deep drawing appliance with the soft side facing up. The foil is conditioned with the bonding agent cloth immediately before heating.

2

2. Bonding agent cloths for conditioning the soft side of the film are separately packed for single-use only.

3

3. The laminated films can now be separated from the model and trimmed to their final length. The finished splints should be about 1-2 mm above the marked equator in order to offer an optimum grip. Tightly fitting splints can easily be shortened – but not lengthened.

Notice for finishing splints

1. The acrylic parts used for finishing are also laminated and do not need to be either finished or polished.
2. The shortening of the deep drawn film and final length of the laminated splints is done with a triangular cutting burr (yellow mark). Smoothing of the machined away, often fibrous splint edges is done with finishing discs.
3. Exposing the head of the anchorage elements from the deep drawing material is done with a polishing disc.

The completed IST®classic^{NEW} appliance


The IST®classic^{NEW} appliance should be cleaned daily with a toothbrush.

The ISTplus® appliance

- 1 Adjustable connecting bridge
- 2 Adjusting screw
- 3 Bridge screw with an Allen key
- 4 Deep drawing splints for the upper and lower jaw
- 5 Plastic anchoring part

- 1 Adjustable connecting bridge
- 2 Adjusting screw
- 3 Bridge screw with an Allen key
- 4 Space saver
- 5 Splint nut
- 6 Plastic anchoring part

Fixation of fasteners in the lower jaw

- First the fasteners are fixed in the lower jaw with kneadable silicone (e.g. Sil-Tech putty/Scheu putty) and aligned
- Check in the articulator!
- In the case of space problems, chamfer the fastening parts in the rear area
- In the correct position, secure the fasteners with acrylate

Fixation of the fastening parts in the upper jaw

- Fixation of the fastening parts with Sil-Tech putty in the upper jaw (in the articulator)
- Fix the fastening parts with acrylate

Dismantling

- Removal of the guiding splints, Allen screws and the space holders


Construction of frontal bite blocks in the lower jaw

1. To protect the fastening parts in the upper jaw, use a strip of insulating foil
2. Build up the bite blocks with acrylate in the lower jaw
3. The bite blocks should be at least as large as the fastening parts in the upper jaw and lie just below it

Finalising

- A1. Do not lift the film from the model. Finalise acrylate on the model.
- A2. The frontal bites are finalised in the lower jaw. Pre-clean the fastening parts in the upper jaw/ lower jaw. Occlusion control in the articulator (see photo. For laminating you need about 2 mm space between the lower jaw bite block and the upper jaw attachment part and more than 2 mm in the molar area).

Lamination

- Embed the upper jaw model in the granulate
- The model should be embedded as deeply as possible (see photo)
- Clamp the LAMItwo film with the soft side up

Lamination of the upper jaw film

- Conditioning of the film with the enclosed adhesive agent cloth. The film is subsequently deep drawn
- Occlusion control in the articulator
- If necessary, grind the bite locks in the lower jaw to make room for the lamination

Lamination of the lower jaw film

- Embed the model as described
- Clamp and condition film
- Deep draw

Finishing the splints*

Lift from the models and finish, expose the screws holes with a spiral drill, mounting of the metal parts.

* The tools required for finishing are offered as a finishing set (Art. no. 98 700).

Final assembly

Activate but how?

In order to ensure a uniform direction of activation of the connecting bridges, it is necessary to observe the following in the case of the final assembly:

- The laser marking must point towards the mesial on the right hand side and towards the distal on the left hand side
- Direction of rotation from the upper jaw to the lower jaw. The connecting bridges shorten (protrusion)
- Direction of rotation from the lower jaw to the upper jaw. The connecting bridges become longer (retrusion)

The ISTplus® appliance should be cleaned daily with a toothbrush.

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