## Das ELASTO-KFO-System

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Addition cross-linking silicone

### Instructions for use

For the production of **SILASTO®** orthodontic appliances

## The ELASTO KFO System

Successful treatment with the ELASTO ORTHODONTIC system requires a silicone-based material with a high level of resilience in case of biologically tenable forces, sustained permanent elasticity, high tensile strength, stable transparency, neutral taste and odour and disinfecting ability. These requirements are met by the silicone material **SILASTO®** 30, 50 and 70 developed by Bayer Dental, Dormagen in collaboration with Dr. Hinz Dental, Herne.

**SILASTO®** which is available in the Shore hardnesses 30, 50 and 70, is a viscous, transparent, addition-curing silicone base material with a platinum catalyst system.

This catalytic system is characterised by high biocompatibility, which makes the material **SILASTO®** physiologically safe to use.

The technical possibilities of a practice laboratory are sufficient for the processing. **SILASTO® 30, 50** and **70** is supplied in 50 ml double chamber cartridges.

#### **SILASTO®**

is mixed on a mixing pot for **SILASTO®** with a mixing gun without a mixing cannula.

The contents of the dual chamber cartridge are sufficient for at least one **SILASTO®** device.

Please observe strict cleanliness when mixing **SILASTO®** and during further processing.

Under no circumstances should **SILASTO®** be allowed to come into contact with tobacco smoke, rubber gloves or organic tin and/or nitrogen compounds when mixing and during further processing. Such contact attacks the catalytic system and affects the hardening of **SILASTO®**.

#### **SILASTO® is further processed using the cuvette technology known from prosthetics as follows:**

- Insulate one cuvette with silicone cream.
- Embed the modelled appliance in the cuvette with the upper jaw pushed down to the wax edge of the lower jaw.
- Spray the cuvette half with the special insulation fluid and produce the counter impression.
- After the plaster setting time, briefly pre-heat the cuvette.
- Open the cuvette carefully and boil out without grease.
- Dry the cuvette halves and allow them to cool down, then spray with the special insulation fluid, allow to work for 30 seconds, remove the excess with an air jet. Then allow to dry for a further 15 minutes.
- Allow the cuvettes to cool down again.
- Fill the cuvette halves with the mixed **SILASTO®**, close the cuvette halves, place these in the screw press and press slowly.
- As soon as material escapes, interrupt the pressing procedure, then continue in intervals in this way until the cuvette edges are tightly closed. The step-wise pressing procedure should not take less than 10 minutes.
- Remove the cuvette from the press and immediately clamp in a cuvette holder.
- Place the cuvette in a pre-heated oven for 45 minutes at 250°C in order to vulcanise the **SILASTO®**.
- Allow the cuvette to cool down and open it carefully.
- Cut off excess pressing material and place the appliance in an ultrasound cleaning bath.
- Process and smooth any remaining unevenness with grinding heads for **SILASTO®** and place the appliance in the ultrasound bath again for a short time.

## SILASTO® polishing fluid

**When using the SILASTO® polishing fluid, please use the special mixing cannula supplied.**

It is important that the mixing cannula remains on the polishing fluid cartridge after use.

Immediately before the next application, the mixing cannula may be replaced with a new one.

Thinly apply the mixed **SILASTO®** polishing fluid with a small, medium-hard brush to the prepared and thoroughly cleaned appliance as follows:

1. Coat from occlusal/incisal appliance with very little polishing fluid (cusps must not taper).
2. Coat the appliance from lingual.
3. Insert clamp from labial into the air slots.
4. Hold the appliance by the clamp and coat it from labial.
5. Hang the appliance thus coated with the clamp on the stand for the second vulcanisation and leave it in a dust-free environment for approx. 5 minutes, so that the brush strokes disappear and a homogeneous, high-gloss surface is created.
6. Place the stand with the appliance hanging on it in the preheated oven and vulcanise at 180 ° C for 10 minutes.

#### **Further materials from our range of Elasto products:**

Order No.	Article
97 459	Polishing fluid mixing cannula (10 pc)
97 114	Special insulation fluid, 250 ml
97 115	Silicone cream, 50 ml
97 216	Filling spatula, straight
97 217	Filling spatula, angled
97 255	Mixing block, DIN A4, 20 sheet
97 256	Mixing block, DIN A5, 20 sheet
97 111	Mixing glass (heat resistant)
97 401	Mixing gun
97 060	Air slot wedges (10 pc)

#### Also available:

Grinding head holders / grinding heads in 4 different sizes / 2 grain sizes

## Processing errors

#### **The appliance is bubbly:**

- It was pressed too quickly
- The cuvette clamp was not tightly closed

#### **The appliance is not fully formed:**

- It was pressed too quickly
- The cuvette was still too warm during the filling process
- Not enough **SILASTO®** material was put in the cuvette

#### **The appliance is sticky:**

- The two material components were not mixed properly
- The vulcanisation time was too short
- The oven temperature was too low
- The **SILASTO®** material processed came into contact with catalyst poisons (e.g., tobacco smoke, hand sweat, hand cream, rubber gloves)

#### **The finished appliance is not transparent:**

- The cuvette was not dry enough before the filling process
- The **SILASTO®** material processed came into contact with water vapor during processing
- The **SILASTO®** material was stored in a damp environment after mixing
- The wrong insulation fluid was used

#### **The appliance has visible seams:**

- The cuvette was still warm during the filling process
- It was pressed too quickly
- The mixed **SILASTO®** material was stored too warm
- The **SILASTO®** material was not processed directly after mixing

#### **There are impurities in the appliance:**

- After the wax was boiled out, processing was not clean
- The wrong insulation fluid was used

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